

BEST AVAILABLE COPY

EXHIBIT A




GL Biochem (Shanghai) Ltd
吉来生化 (上海) 有限公司

Product Specification

Name: **DMT-Cl**

Category: N-Protecting Reagents



Product Data Sheet

Product Name	DMT-Cl; 4,4'-Dimethoxytrityl chloride; 4,4'-Dimethoxytriphenylmethyl chloride
CAS No.	40615-36-9
Molecular Formula	C ₂₁ H ₁₉ ClO ₂
Molecular Weight	338.8
Appearance	Light pink color powder
Purity (HPLC)	98% min.
Melting Point	120-125 °C (dec.)
TLC Analysis	One spot
300 MHz ¹ H NMR Spectrum (CDCl ₃)	Consistent
Loss on Drying	< 0.5%
Solubility Test (In Pyridine)	Clear solution with 1g/10ml
Solvent of Recrystallization	Benzene, Hexanes
Use Test	Passed

EXHIBIT B

PROTECTIVE GROUPS IN ORGANIC SYNTHESIS

Second Edition

THEODORA W. GREENE

The Rowland Institute for Science, Inc.

and

PETER G. M. WUTS

The Upjohn Company



A WILEY-INTERSCIENCE PUBLICATION

JOHN WILEY & SONS, INC.

New York / Chichester / Brisbane / Toronto / Singapore

PREFACE TO THE SECOND EDITION

In recognition of the importance of preserving what has been written, it is a policy of John Wiley & Sons, Inc. to have books of enduring value published in the United States printed on acid-free paper, and we extend our best efforts to that end.

Copyright © 1991 by John Wiley & Sons, Inc.
All rights reserved. Published simultaneously in Canada.

Reproduction or translation of any part of this work beyond that permitted by Section 107 or 108 of the 1976 United States Copyright Act without the permission of the copyright owner is unlawful. Requests for permission or further information should be addressed to the Permissions Department, John Wiley & Sons, Inc.

Library of Congress Cataloging in Publication Data:

Greene, Theodore W., 1931-
Protective groups in organic synthesis / Theodore W. Greene and Peter G. M. Wuts.—2nd ed.
p. cm.

Includes index.

ISBN 0-471-62304-6

1. Organic compounds—Synthesis. 2. Protective groups (Chemistry)

I. Wuts, Peter G. M.

QA261.G665 1991

547.2—dc20

90-45704
CIP

Printed in the United States of America

10 9 8 7 6

Since publication of the first edition of this book in 1981, many new protective groups and many new methods of introduction or removal of known protective groups have been developed: 206 new groups and approximately 1500 new references have been added. Most of the information from the first edition has been retained. To conserve space, generic structures used to describe Formation/Cleavage reactions have been replaced by a single line of conditions, sometimes with explanatory comments, especially regarding selectivity. Some of the new information has been obtained from Online Searches of Chemical Abstracts, which have limitations. For example, Chemical Abstracts indexes a review article about protective groups only if that word appears in the title of the article. References are complete through 1989. Some references, from more widely circulating journals, are included for 1990.


Two new sections on the protection for indoles, imidazoles, and pyrroles, and protection for the amide —NH are included. They are separated from the regular amines because their chemical properties are sufficiently different to affect the chemistry of protection and deprotection. The Reactivity Charts in Chapter 8 are identical to those in the first edition. The chart number appears beside the name of each protective group when it is first discussed.

A number of people must be thanked for their contributions and help in completing this project. I am grateful to Gordon Bundy, who loaned me this card file, which provided many references that the computer failed to find, and to Bob Williams, Spencer Knapp, and Tadahiko Fukuyama for many references on amine and amide protection. I thank Theo Greene who checked and rechecked the manuscript for spelling and consistency and for the herculean task of checking all the references to make sure my 3s and 8s and 7s and 9s were not interchanged, all without

12. Et_2AlCl , CH_2Cl_2 , 3 min, 70–85% yield.¹⁹ This method was used to remove the trityl group from various protected deoxyribonucleosides. The TBDPS group is stable to these conditions.

13. TsOH , MeOH , 25° , 5 h.²⁰

14. $\text{BF}_3 \cdot \text{Et}_2\text{O}$, $\text{HSCH}_2\text{CH}_2\text{SH}$, 80% yield.²¹



15. Na , NH_3 ,²² Benzyl groups are also removed under these conditions.

16. ZnBr_2 , MeOH , 100% yield.²³

57. α -Naphthyl(diphenylmethyl) Ether: $\text{RO}-\text{C}(\text{Ph})_2-\alpha\text{-C}_{10}\text{H}_7$ (Chart 1)

The α -naphthyl(diphenylmethyl) ether was prepared to protect, selectively, the 5'-OH group in nucleosides. It is prepared from α -naphthyl(diphenylmethyl) chloride in pyridine (65% yield), and cleaved selectively in the presence of a *p*-methoxyphenyldiphenylmethyl ether with sodium anthracene, a (THF, 97% yield). The *p*-methoxyphenyldiphenylmethyl ether can be cleaved with acid in the presence of this group.²⁴



58. *p*-Methoxyphenyldiphenylmethyl Ether (MMT-OR):
 $p\text{-MeOC}_6\text{H}_4(\text{Ph})_2\text{C-OR}$ (Chart 1)

59. Di(*p*-methoxyphenyl)phenylmethyl Ether (DMT-OR):
 $(p\text{-MeOC}_6\text{H}_4)_2\text{PhC-OR}$

60. Tri(*p*-methoxyphenyl)methyl Ether (TMT-OR): $(p\text{-MeOC}_6\text{H}_4)_3\text{C-OR}$

These were originally prepared by Khosla²⁵ as selective protective groups for the 5'-OH of nucleosides and nucleotides. They were designed to be more acid-labile than the trityl group because depurination is often a problem in the acid-catalyzed removal of the trityl group. Introduction of *p*-methoxy groups increases the rate of hydrolysis by about one order of magnitude for each *p*-methoxy substituent. For 5'-protected uridine derivatives in 80% AcOH , 20° , the time for hydrolysis was

as follows:

$(p\text{-MeOC}_6\text{H}_4)_n(\text{Ph})_m\text{C-OR}$

$n = 0, m = 3, 48 \text{ h}$

$n = 1, m = 2, 2 \text{ h}$

$n = 2, m = 1, 15 \text{ min}$

$n = 3, m = 0, 1 \text{ min}$

The trimethoxy derivative is too labile for most applications, but the mono and di-derivatives have been used extensively in the preparation of oligonucleotides and oligonucleosides. The monomethoxy derivative has been used for the selective protection of a primary allylic alcohol over a secondary allylic alcohol (MMT, Pyr, -10°).²⁶

Cleavage


In practice the various trityl derivatives are cleaved with acid, but the monomethoxy derivative can be cleaved with sodium anthracene in HMFA (90% yield).²⁷ It is not cleaved by sodium anthracene, used to cleave α -naphthyl(diphenylmethyl) ethers.²⁴

A solution of 3% $\text{CCl}_3\text{CO}_2\text{H}$ in 95:5 $\text{CH}_3\text{NO}_2/\text{MeOH}$ is recommended for removal of the DMT group from the 5'-OH of deoxyribonucleosides because of reduced levels of depurination compared to $\text{Cl}_3\text{CO}_2\text{H}/\text{CH}_2\text{Cl}_2$, $\text{PhSO}_3\text{H}/\text{MeOH}/\text{CH}_2\text{Cl}_2$, and $\text{ZnBr}_2/\text{CH}_3\text{NO}_2$.²⁸

61. 4-(4'-Bromophenoxyphenyl)diphenylmethyl Ether:
 $p\text{-(}p\text{-BrC}_6\text{H}_4\text{C(OCH}_2\text{C}_6\text{H}_4)_2\text{C}_6\text{H}_4\text{)(Ph)}_2\text{C-OR}$

This group was developed for protection of the 5'-OH group in nucleosides. The derivative is prepared from the corresponding triarylmethyl chloride, and is cleaved by reductive cleavage (Zn/AsOH) of the phenacyl ether to the *p*-hydroxyphenyldiphenylmethyl ether followed by acidic hydrolysis with formic acid.²⁹

62. 4,4'-Tribis(4,5-dichlorophthalimido)phenylmethyl Ether (CPT-OR):



The CPT group was developed for the protection of the 5'-OH of ribonucleosides. It is introduced with $\text{CPTNa}/\text{AgNO}_3/\text{DMF}$ (15 min) in 80–96% yield and can be removed by ammonia followed by 0.01 *M* HCl or 80% AcOH .³⁰ It can also be removed with hydrazine and acetic acid.³¹

**This Page is Inserted by IFW Indexing and Scanning
Operations and is not part of the Official Record**

BEST AVAILABLE IMAGES

Defective images within this document are accurate representations of the original documents submitted by the applicant.

Defects in the images include but are not limited to the items checked:

☒ **BLACK BORDERS**

☐ **IMAGE CUT OFF AT TOP, BOTTOM OR SIDES**

☒ **FADED TEXT OR DRAWING**

☐ **BLURRED OR ILLEGIBLE TEXT OR DRAWING**

☐ **SKEWED/SLANTED IMAGES**

☐ **COLOR OR BLACK AND WHITE PHOTOGRAPHS**

☐ **GRAY SCALE DOCUMENTS**

☐ **LINES OR MARKS ON ORIGINAL DOCUMENT**

☐ **REFERENCE(S) OR EXHIBIT(S) SUBMITTED ARE POOR QUALITY**

☐ **OTHER:** _____

IMAGES ARE BEST AVAILABLE COPY.

As rescanning these documents will not correct the image problems checked, please do not report these problems to the IFW Image Problem Mailbox.